

# Cerebrovascular Reserve Testing – A case with adequate reserve

## PHYSICIAN'S CONCERNS?


An active, 55 year old man presents with a left internal carotid artery occlusion and recurring symptoms. He has been referred by his primary care physician.

## HOW DID NOVA HELP?

Quantitative flow measurements pre and post diamox are used to evaluate cerebrovascular reserve and rule out cerebral ischemia as the origin of symptoms.

## Patient History

- ❖ 55 year old right handed man
- ❖ Left internal carotid artery occlusion, likely dissecting in origin
- ❖ Initial presentation 13 months ago when he had an episode of left facial hypersensitivity, headaches and droopy eyelid; all consistent with Horner's syndrome.
- ❖ Past medical history relevant for hypertension and hyperlipidemia
- ❖ Family history of coronary artery disease (father died at 49)
- ❖ High level of cognitive function; employed as an engineer
- ❖ Residual Horner's syndrome. Currently, he is not symptomatic for TIA or stroke



**Figure 1:** NOVA 3D with slice plane on left middle cerebral artery

## NOVA Report

- Flow measurements in the major cerebral vessels were taken at baseline and following administration of the vasodilator, diamox.
- Blood flow increased in response to diamox in all measured vessels.
- The flows in the L and R MCAs increased by 36% and 52% respectively, demonstrating good cerebrovascular reserve. (Figure 2)

# Cerebrovascular Reserve Testing – A case with adequate reserve




Figure 2: NOVA Vessel Map

## Conclusions

- NOVA qMRA quantifies the collateral flow supply
- In response to a diamox challenge, the patient demonstrates good cerebrovascular reserve
- Patient is adequately compensated and does not need flow augmentation